

# New. Perfect.

The Degidi Plugger from Ustomed

[www.ustomed.de](http://www.ustomed.de)



Biomaterial plugger  
for socket preservation

 SLINE




Competence in Instruments  
MADE IN GERMANY

## Degidi Plugger Precise. Perfect.

In case of immediate implantation, it is very important to fill the gap between the implant and the buccal plate with some low resorbable biomaterial.

In order to pack this space in the best possible way we have designed 2 instruments that due to the different shapes (4 different working ends) and dimensions of the apical wedge, fit perfectly in this gap.



- 1 | Placed Implant.
- 2 | Application of biomaterial.
- 3 | Condensing of biomaterial with 4 different sizes.
- 4 | Result.

## Degidi Plugger Functionality. Design.

The Degidi Plugger from USTOMED is an ergonomically designed instrument for the socket preservation.

4 different and very accurate shaped working ends allow the surgeon to place, fill and condense the biomaterial always in the perfect way.





### Dr. Marco Degidi

Adjunct Professor, Postgraduate Program in Clinic Implantology and Biomaterials, University of Chieti. Adjunct Professor, Postgraduate Master in Clinic Implantology, University of Bologna.

Speaker in many national and international seminars, congresses and courses; many national and international publications.

Active member of AO, EAO, ICOI and SIO.


Private Practice in Bologna, restricted to implant surgery & prosthetic.



# Degidi Plugger

Size 2

Size 4



1/1

Order Number  
67-693-023

1/1

Order Number  
67-693-014

Size 3

Size 1

 **USTOMED**<sup>®</sup>  
INSTRUMENTE

USTOMED INSTRUMENTE  
Ulrich Storz GmbH & Co. KG

Am Eichbuehl 70 – 72  
78532 Tuttlingen / Germany

Hersteller von Chirurgie- und  
Dentalinstrumenten

Phone: +49 (0) 74 61 / 9 65 85 0  
Fax: +49 (0) 74 61 / 9 65 85 65

Manufacturer of Surgical and  
Dental Instruments

info@ustomed.de  
www.ustomed.de